



[UserGuide](#) /

CropDicom


User Guide: CropDicomTool

last update: 2011/4/20


- The idea of the CropDicomTool is to prepare the images for the harris lines tool (HL tool)
- Each input image for the HL tool has to be in the dicom format and contain one single bone, which is vertically oriented
- the buttons on the right hand side are used from top to down

(1) open dicom or tiff image

- black background needed, otherwise 'invert' the image
- currently only tiff and dicom images can be loaded (other formats upon request [Contact](#))




(1) open file



opening panel

(2) extract bones

- threshold the input image such that you see your bones in white on the right image
- calculate the automatic threshold (with button)
- tune / adjust the calculated threshold with the offset (play with negative or positive small offsets)
- every bone should be clearly separated from other things, i.e., no connecting white line should be visible between the bones
- small numbers or box artefacts should also be separated, they can be removed later



(2) thresholded bones

(2) extract bones

compute automatic threshold

threshold


offset

apply threshold

threshold panel

(3) separate bones

- one color per bone
- mark each single bone with a separate color (play with the input field to increase or decrease the number)
- e.g. by increasing the value (min area size), small numbers or similar things can be removed



(3) labeled bones (one color per bone)

(3) separate bones

enter min area size

automatic labeling /coloring

labeling panel

(4) save single bones

- one dcm-file per bone
- each colored bone is extracted and saved
- you have to provide a file name for every bone
- before saving, each image is displayed such that you can give the correct name to it
- artifacts are also saved, but they can be manually deleted after the saving procedure
- the bone should be oriented vertically (if this is not the case, it needs to be rotated in another tool, rotation is currently not provided in the CropDicomTool)